wiroCOMPLETE KNIT

| | | |
|------|--|-------------|
| [de] | Produktdatenblatt PP+PE KITTEL | Seite 2 |
| [en] | Product data sheet PP+PE GOWN | Page 6 |
| [nl] | Productgegevensformulier PP+PE JAS | Zijde 10 |
| [fr] | Fiche produit BLOUSES PP+PE | Page 14 |
| [es] | Hoja de datos del producto BATA DE PP+PE | Página 18 |
| [it] | Scheda prodotto CAMICE IN PP+PE | Pagina 22 |
| [el] | Δελτίο δεδομένων προϊόντος ΙΑΤΡΙΚΗ ΠΟΔΙΑ PP+PE | Σελίδα 26 |
| [tr] | Ürün Bilgi Formu PP+PE ÖNLÜK | Sayfa 30 |
| [ru] | Техническая спецификация ХАЛАТ ИЗ ПП+ПЭ | Сторона 34 |
| [pl] | Karta produktu KITEL PP+PE | Strona 38 |
| [da] | Produktdatablad PP + PE-KITTEL | Side 42 |
| [sv] | Produktdatablad PP+PE-SKYDDSROCK | Sida 46 |
| [no] | Produktdatablad PP+PE-KITTEL | Side 50 |
| [cs] | Informační list výrobku PP+PE PLÁŠŤ | Strana 54 |
| [hu] | Termékadatlap PP+PE LÁTOGATÓI KÖPENY | Oldal: 58 |
| [sk] | Informačný list výrobku PP + PE PLÁŠŤ | Strana 62 |
| [sl] | Podatkovni list izdelka PP+PE HALJA | Stran 66 |
| [ro] | Fișa produsului HALAT PP+PE | Pagina 70 |
| [pt] | Folha de dados do produto BATA PP+PE | Página 74 |
| [lt] | Gaminio duomenų lapas PP+PE DARBINIAI CHALATAI | Puslapis 78 |
| [lv] | Produkta datu lapa PP+PE UZSVĀRČĪ | Lappuse 82 |
| [bg] | Информационен лист за продукта PP+PE ПРЕСТИЛКА | Страница 86 |
| [uk] | Технічний лист ХАЛАТ З ПП+ПЕ | Сторінка 90 |
| [et] | Toote teabeleht PP+PE KITTEL | Lk 94 |
| [fi] | Tuotetietolehti PP + PE -TAKIT | Sivu 98 |

wiroCOMPLETE KNIT

PP+PE KITTEL

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wiroCOMPLETE KNIT

PP+PE GOWN

PE full coating | Knitted cuffs | Neck & waist ties

PROFESSIONAL

Product data sheet

wiroCOMPLETE KNIT

Manufacturer

| | |
|-------------------|---|
| Company + Address | WIROS Wilfried Rosbach GmbH, Karl-Arnold-Str. 5, 47877 Willich, Germany |
| SRN | DE-MF-000004956 |
| Trademark | care & serve® |

Authorisation

| | |
|--------|--|
| CH REP | Swiss AR Services AG, Industriestrasse 47, 6300 Zug, Switzerland |
|--------|--|

Product description


| | |
|---------------------|--|
| Category | PROFESSIONAL |
| Product description | care & serve® PROFESSIONAL wiroCOMPLETE KNIT PP+PE gown PE full coating Knitted cuffs Neck & waist ties |
| Product features | PE full coating Knitted cuffs Neck & waist ties |

Material

| | |
|-----------------|--|
| Material | Polypropylene spun-bonded fabric with polyethylene coating |
| Material weight | PP: $\approx 18 \text{ g/m}^2 (\pm 2 \text{ g/m}^2)$ |

Application

| | |
|---------------------------|--|
| Intended use | Resistant medical protective clothing as defined in the Technical Rules for Biological Agents (TRBA 250) to prevent contamination of work or normal clothes during interventions or care measures resulting from direct contact with blood, secretions, excretions or other contaminated material. The material is tested for resistance to penetration by contaminated liquids, infectious agents, contaminated solid articles and contaminated liquid aerosols in accordance with 4.1.4 EN 14126:2003 + AC:2004. |
| Contact with food produce | Short-term indirect contact with foodstuffs is permitted. |
| Scope of application | There is no direct medical application on the patient. Medical protective clothing is intended to cover clothing and/or hair. |
| Application group | Employees and visitors from laboratories, the pharmaceutical industry, hospitals or similar institutions. |
| Patient group | Medical protective clothing is suitable for all patient groups. |
| Indication | Medical clothing as a barrier to cover clothing and/or hair. |
| Contraindication | No known contraindications. |
| Notes on use | Disposable product Not sterile |
| Limitation | No PPE (personal protective equipment) Not suitable for packing or storing food. |
| Note | All serious incidents directly related to this product must be reported to the manufacturer and the competent authorities. |
| Application areas | Hospital Old people's home inc. specialist circles Medical practices Industry Pharmaceutical industry Laboratory Food processing Cleaning Hygiene |

Product data sheet

wiroCOMPLETE KNIT

| | | |
|-----------------------------|--|---|
| Supplementary material info | Latex-free |  |
| Health advice | If skin irritation or allergic reactions occur, discontinue use immediately and seek medical advice. As far as is known according to the current state of the art, the products do not contain any toxic, carcinogenic, mutagenic or reproductively toxic substances when used in accordance with their intended use. These products do not contain latex components and are considered hypoallergenic. However, they can cause irritation in sensitive skin. Long-term wear may restrict skin ventilation, which can cause skin irritation. | |

Standards

| CE marking | Class I medical device in accordance with EU 2017/745 |   | | | | | | | | | | |
|--|---|---|-------|--|---|--|---|---|---|--|---|--|
| SRN | DE-MF-000004956 | | | | | | | | | | | |
| Basis UDI DI | 4051642-wiroCOMPLETEKN-5U | | | | | | | | | | | |
| | CE marking The CE mark guarantees the free movement of products and goods within the European Economic Community. All products with the CE mark comply with the essential requirements of the European regulations specified. | | | | | | | | | | | |
| EN 14126:2003 - 4.1.4 | 4.1.4 Performance requirements for resistance to infectious agent penetration | | | | | | | | | | | |
| | <table><thead><tr><th>Resistance to penetration</th><th>Class</th></tr></thead><tbody><tr><td>Contaminated fluids according to ISO 16603 and ISO 16604</td><td>4</td></tr><tr><td>From infectious agents according to EN ISO 22610</td><td>6</td></tr><tr><td>Contaminated liquid aerosols according to ISO 22611</td><td>1</td></tr><tr><td>Contaminated solid particles according to EN ISO 22612</td><td>3</td></tr></tbody></table> | Resistance to penetration | Class | Contaminated fluids according to ISO 16603 and ISO 16604 | 4 | From infectious agents according to EN ISO 22610 | 6 | Contaminated liquid aerosols according to ISO 22611 | 1 | Contaminated solid particles according to EN ISO 22612 | 3 | |
| Resistance to penetration | Class | | | | | | | | | | | |
| Contaminated fluids according to ISO 16603 and ISO 16604 | 4 | | | | | | | | | | | |
| From infectious agents according to EN ISO 22610 | 6 | | | | | | | | | | | |
| Contaminated liquid aerosols according to ISO 22611 | 1 | | | | | | | | | | | |
| Contaminated solid particles according to EN ISO 22612 | 3 | | | | | | | | | | | |
| | Information The tests according to EN 14126:2003 + AC:2004 Point 4.1.4 only comprise one part of the requirements of the whole EN 14126:2003 + AC:2004. Other requirements were neither checked nor confirmed, it is therefore not infectious agent protective clothing according to (EU) 2016/425 PPE and must therefore not be labelled as such. | | | | | | | | | | | |
| VO (EU) 2023/988 | Regulation (EU) No 2023/988 on general product safety | | | | | | | | | | | |
| VO (EG) 1935/2004 | Regulation regarding plastic materials and articles intended to come into contact with food. |  | | | | | | | | | | |
| VO (EU) 10/2011 | Regulation on plastic materials and articles intended to come into contact with food | | | | | | | | | | | |

Use

| | | |
|----------------------|--|---|
| Instructions for use | Before use Please make sure that you have selected the right product for your needs and areas of use. Please do not wear medical protective clothing if it is faulty or has already been used. Medical protective clothing is intended for short-term use and use is limited to one patient in the medical field. Use Gowns are worn with the closed side at the front and are fastened at the back or side. To do this, slide your arms into the sleeves and fasten the ties on your neck and waist. In order to prevent contamination with clothing when taking the gown off, it is recommended for a second person to untie the gown from behind and slide it forwards over your shoulders. The person wearing the gown grasps the upper ends of the gown and turns it inside out when removing it. | |
| Safety notice | According to our risk analysis under EN ISO 14971:2019 + A11:2021, no additional safety instructions are required for this product. | |
| Disposal | Unused and uncontaminated products can be incinerated in an environmentally friendly manner or disposed of at a waste disposal site. Contaminated products must be disposed of in accordance with the national laws and regulations valid for the respective type of contamination. | |
| Shelf life | 5 years from the production date, if unused and stored correctly. | |
| Storage | Keep out of direct sunlight. Store in a dry place. Incorrect storage can lead to a reduction in the shelf life. |   |

Product data sheet

wiroCOMPLETE KNIT

Traceability Batch-related traceability from the determination of requirements to issuing.
Product data sheet Version 2 Valid from Batch 2617 +

Packaging

Recycling Licensed in accordance with VerpackG (German Packaging Act), Licence number DE1618662887322.
Please recycle the packaging.

www.wiros.de

WEB <https://www.wiros.de/shop/kittel/pppe-kittel/wirocomplete-knit/>

Variant 1

Sales unit 100 Unit/Cardboard box
Packaging level 1 10 Unit/Polybag
Minimum order quantity 1 Cardboard box

| Item numbers | Item number | Colour | Size | Size (circumference x length) | Tolerance |
|--------------|-------------|--------|------|-------------------------------|-----------|
| | 25305744 | yellow | XL | ≈ 140 x 135 cm | ± 3 cm |

| GTINs | Item number | Unit of Use UDI | Packaging level 1 UDI | Packaging level 2 UDI | Cardboard box UDI |
|-------|-------------|--------------------|--------------------------|--------------------------|----------------------|
| | 25305744 | 4051642004877 | 4051642075419 | ----- | 4051642575414 |

The numerical data for the measurements, weights, thicknesses and tolerances are target values, and there may be some deviation.
The information represents the current state of our knowledge. It offers no guarantee of the product characteristics and does not form the basis of a contractual legal understanding. Please reconfirm any information that is older than 2 years.
Typographical and other errors excepted.
The product data sheets provided are the property of WIROS Wilfried Rosbach GmbH. WIROS Wilfried Rosbach GmbH is the owner of all copyrights and rights of use and exploitation for the product data sheets. It is not permitted to change the product data sheets or to use the pictograms. The product data sheets provided may only be used in conjunction with original WIROS products.

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