wiroMICROGRIP PF

[de]	Produktdatenblatt LATEX HANDSCHUHE	Seite 2
[en]	Product data sheet LATEX GLOVES	Page 8
[nl]	Productgegevensformulier LATEX HANDSCHOENEN	Zijde 14
[fr]	Fiche produit GANTS EN LATEX	Page 20
[es]	Hoja de datos del producto GUANTES DE LÁTEX	Página 26
[it]	Scheda prodotto GUANTI IN LATTICE	Pagina 32
[el]	Δελτίο δεδομένων προϊόντος ΓΑΝΤΙΑ ΛΑΤΕΞ	Σελίδα 38
[tr]	Ürün Bilgi Formu LATEKS ELDIVENLER	Sayfa 44
[ru]	Техническая спецификация ПЕРЧАТКИ ЛАТЕКСНЫЕ	Страница 49
[pl]	Karta produktu RĘKAWICE LATEKSOWE	Strona 55
[da]	Produktdatablad LATEX-HANDSKER	Side 60
[sv]	Produktdatablad LATEXHANDSKAR	Sida 65
[no]	Produktdatablad LATEX-HANSKER	Side 70
[cs]	Informační list výrobku LATEXOVÉ RUKAVICE	Strana 75
[hu]	Termékadatlap LATEX VÉDŐKESZTYŰ	Oldal: 80
[sk]	Informačný list výrobku LATEXOVÉ RUKAVICE	Strana 85
[sl]	Podatkovni list izdelka ROKAVICE IZ LATEKSA	Stran 90
[ro]	Fișa produsului MÂNUȘI DIN LATEX	Pagina 95
[pt]	Folha de dados do produto LUVAS DE LÁTEX	Página 100
[lt]	Gaminio duomenų lapas LATEKSO PIRŠTINĖS	Puslapis 105
[lv]	Produkta datu lapa LATEKSA CIMDI	Lappuse 110
[bg]	Информационен лист за продукта ЛАТЕКСОВИ РЪКАВИЦИ	Страница 115
[uk]	Технічний лист ЛАТЕКСНІ РУКАВИЧКИ	Сторінка 120
[et]	Toote teabeleht LATEXSKINDAD	Lk 125
[fi]	Tuotetietolehti LATEKSIKÄSINEET	Sivu 130

LATEX HANDSCHUHE

Puderfrei | Micro-texturiert

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LATEX GLOVES

Powder free | Micro-textured

i* This product contains natural rubber latex, which can cause allergic reactions including anaphylactic reactions.

PREMIUM

Manufacturer

Company + Address	WIROS Wilfried Rosbach GmbH, Karl-Arnold-Str. 5, 47877 Willich, Germany
SRN	DE-MF-000004956
Trademark	care & serve®

Authorisation

CH REP	Swiss AR Services AG, Industriestrasse 47, 6300 Zug, Switzerland
--------	--

Product description

Category	PREMIUM
Product description	care & serve® PREMIUM wirosMICROGRIP PF Latex gloves Powder free Micro-textured i* This product contains natural rubber latex, which can cause allergic reactions including anaphylactic reactions.

Product features	Powder free Micro-textured
------------------	-------------------------------

Material

Material	Natural latex
Thickness	Finger: $\approx 0,13$ mm ($\pm 0,01$ mm) Palm: $\approx 0,10$ mm ($\pm 0,01$ mm) Hand edge: $\approx 0,08$ mm ($\pm 0,01$ mm)
Material weight	Size M $\approx 5,8$ g/pc. ($\pm 0,2$ g)

Application

Intended use	Disposable protective glove to avoid cross-contamination between users and patients when performing medical examinations, and to prevent the surface of the hands from coming into contact with substances, mixtures, or harmful biological agents that are hazardous to health. The level of protection and the maximum wearing time come from the determined performance level for each material tested.
Contact with food produce	Short-term indirect contact with foodstuffs is permitted.
Scope of application	Gloves are used to prevent contamination of the hands of medical staff.
Application group	Use is restricted to appropriately trained professionals. Not suitable for users with a latex allergy.
Patient group	Not suitable for patients with a latex allergy. The use of Latex gloves in neonatal care is not recommended in order to prevent sensitisation.
Indication	Gloves are used to prevent contamination of the hands of medical personnel, 1) in case of direct contact with blood or body fluids, mucous membranes, non-intact skin and other potentially infectious material. 2) in direct contact with patients colonised or infected with pathogens transmitted by the contact route, e.g. VRE, MRSA, RSV; or 3) when handling or touching visibly or potentially contaminated patient care equipment and environmental surfaces.
Contraindication	Latex can cause allergic reactions and sensitisation.
Notes on use	Disposable product Not sterile

Limitation	It is the sole responsibility of the user to determine the suitability and/or duration of use of care & serve® products for their intended use. The manufacturer assumes no responsibility for the improper use of care & serve® products. If you have a known allergy to latex, please do not use this product. If a risk analysis has deemed that the protection listed above is NOT sufficient, protective clothing with higher protective properties should be used. Not suitable for use with substances which contain petroleum or oil. Not suitable for packing or storing food.		
Note	All serious incidents directly related to this product must be reported to the manufacturer and the competent authorities.		
Application areas	Hospital Old people's home inc. specialist circles Medical practices Industry Pharmaceutical industry Laboratory Food processing HACCP Service Food & drink Hospitality Cleaning Hygiene Cosmetics Wellness Beauty industry House + Garden Agriculture Construction Industry		
Health advice	Contains natural rubber latex that can cause type IV and type I allergies. If skin irritation or allergic reactions occur, discontinue use immediately and seek medical advice. As far as is known according to the current state of the art, the products do not contain any toxic, carcinogenic, mutagenic or reproductively toxic substances when used in accordance with their intended use.		

Standards

CE marking

SRN

Basis UDI DI

Class I medical device in accordance with EU 2017/745

Personal protective equipment category III in accordance with Regulation (EU) 2016/425

DE-MF-000004956

4051642-wiroMICROGRIPP-7N

CE marking

The CE mark guarantees the free movement of products and goods within the European Economic Community. All products with the CE mark comply with the essential requirements of the European regulations specified.

CE 2777

MD CLASS I

PPE CAT. III

EN ISO 21420:2020

Standard: Protective gloves - General requirements and test methods

Test	Result
pH value	Passed
Mobility	5
Minimum length	Length
Size 7: 230 mm	≈ 240 mm
Size 8: 240 mm	≈ 240 mm
Size 9: 250 mm	≈ 240 mm
Size 10: 260 mm	≈ 240 mm

Information

Observe the product information and instructions for use.

www.wiros.de/IFI

EN 455-1:2022

Medical gloves for single use - Part 1: Requirements and testing for freedom from holes

Water retention test	Result
Seal tightness	AQL 1.5

AQL 1.5

EN 455-2:2015

Medical gloves for single use - Part 2: Requirements and testing for physical properties

Minimum length	Length
XS - XL ≈ 240 mm	≈ 240 mm
Test	Result
Tensile strength	Passed

EN 455-3:2015

Medical gloves for single use - Part 3: Requirements and testing for biological evaluation

Test	Result
Powder residues	< 2 mg
Extractable proteins	< 50 µg/g

Information

i* This product contains natural rubber latex, which can cause allergic reactions including anaphylactic reactions.

Clear labelling of whether the glove is powdered or powder-free is mandatory.

EN 455-4:2009

Medical gloves for single use - Part 4: Requirements and testing for shelf life determination

Test	Result
Shelf life	3

EN ISO 374-1:2016 +
A1:2018 TYPE B

Standard: Protective gloves against dangerous chemicals and micro-organisms - Part 1: Terminology and performance requirements for chemical risks.

Test chemical	Level
K Sodium Hydroxide 40%	4
P Hydrogen Peroxide 30 %	2
T Formaldehyde 37%	3

Information

Breakthrough time: Level 1 > 10 min. | Level 2 > 30 min. | Level 3 > 60 min. | Level 4 > 120 min. | Level 5 > 240 min. | Level 6 > 480 min.

This information does not specify the actual protection duration in the workplace and does not distinguish between mixtures and pure chemicals.

The resistance to chemicals is assessed under laboratory conditions using samples only taken from the palm of the hand (except where the glove measures 400 mm or longer - in this case, the cuff is also tested), and applies only to the tested chemicals. The resistance can alter when the chemicals are used in a mixture.

We recommend checking whether the gloves are suitable for the intended application, as the conditions at the workplace can vary from those in the type test depending on the temperature, abrasion and degradation.

If protective gloves are already in use, these can offer reduced resistance to dangerous chemicals due to changes in their physical properties. Degradation caused by contact with chemicals, movements pulling threads, rubbing, etc. can considerably reduce the actual duration of use. Degradation can be the most important factor with aggressive chemicals, which must be taken into consideration when selecting chemical-resistant gloves.

The gloves must be checked for flaws or defects before use.

For single use only, do not reuse.

EN ISO 374-2:2019

Protective gloves against dangerous chemicals and microorganisms - Part 2: Determination of resistance to penetration

Test	Result
Air penetration test	Passed

EN ISO 374-4:2019

Protective gloves against dangerous chemicals and microorganisms - Part 4: Determination of resistance to degradation by chemicals

Test chemical	Degradation %
K Sodium Hydroxide 40%	-32,8
P Hydrogen Peroxide 30 %	-6,5
T Formaldehyde 37%	-13,7

EN ISO 374-5:2016 VIRUS

Standard: Protective gloves against dangerous chemicals and micro-organisms - Part 5: Terminology and performance requirements for micro-organism risks

Test	Result
Protection against bacteria and fungi	Passed
Protection against bacteria, fungi and viruses	Passed

Information

The penetration resistance is assessed under laboratory conditions and refers solely to the tested samples.

VO (EU) 2023/988

Regulation (EU) No 2023/988 on general product safety

VO (EG) 1935/2004

Regulation regarding plastic materials and articles intended to come into contact with food.

VO (EU) 10/2011

Regulation on plastic materials and articles intended to come into contact with food

Use

Instructions for use

Before use

Please note that this product contains latex. Please make sure that you have selected the correct product for your requirements and application areas, and observe the requirements set out in the respective standards. Should you have any questions, please consult the supplier of your disposable protective gloves. Please do not wear the disposable protective gloves if they show visible defects or have already been used. Disposable protective gloves are only intended for short-term use when working with light chemicals.

Use

Please select the correct size for you.

The disposable protective gloves protect the area they cover. If protection is required for a different part of the body or the respiratory system, use additional PPE products and, for each application, only use the PPE approved for that intended purpose. If medical protection is required for a different part of the body or the respiratory system, use additional medical products and, for each application, only use the medical products approved for that intended purpose. The user is solely responsible for the correct combination of additional PPE or medical products and compliance with the necessary requirements of the standards for each application. You alone are responsible for the duration of use during the respective application, while taking into consideration the protective effect and wear comfort. In the medical sector, use is limited to one patient. An alteration in the given protection level cannot be excluded following contamination. Also check regularly during use that the item is intact. Replace a damaged disposable protective glove with an intact one immediately. Remove the disposable protective glove if it is contaminated or damaged and dispose of it correctly.

Warning: In the event of a potential allergic reaction, stop using this product immediately. Wash your hands thoroughly and seek medical help if necessary.

Safety notice

Contains natural latex that can cause allergic reactions.

Disposal

Unused and uncontaminated products can be incinerated in an environmentally friendly manner or disposed of at a waste disposal site. Contaminated products must be disposed of in accordance with the national laws and regulations valid for the respective type of contamination.

Shelf life

3 years from the production date, if unused and stored correctly.

Storage

Protect from intense heat, sunlight, frost, moisture, ozone, X-rays and ultraviolet light.

Store in a dry place.

Storage temperature between 10 °C and 40 °C.

Incorrect storage can lead to a reduction in the shelf life.

Traceability

Batch-related traceability from the determination of requirements to issuing.

Product data sheet Version 3 Valid from Batch 2617 +

Packaging

Recycling

Licensed in accordance with VerpackG (German Packaging Act), Licence number DE1618662887322.

Please recycle the packaging.

www.wiros.de

WEB

<https://www.wiros.de/shop/handschuhe/latex-handschuhe/wiromicrogrip-pf/>

Variant 1

Sales unit

1000 Unit/Cardboard box

Packaging level 1

100 Unit/Dispenser box

Minimum order quantity

1 Cardboard box

Item numbers

Item number	Colour	Size	Size (W x L)	Tolerance
40351401	white	S	≈ 80 x 240 mm	± 5 mm
40351402	white	M	≈ 90 x 240 mm	± 5 mm
40351403	white	L	≈ 100 x 240 mm	± 5 mm
40351404	white	XL	≈ 110 x 240 mm	± 5 mm

wiroMICROGRIP PF

GTINs	Item number	Unit of Use UDI	Packaging level 1 UDI	Packaging level 2 UDI	Cardboard box UDI
	40351401	4051642005966	4051642150239	-----	4051642650234
	40351402	4051642005973	4051642150246	-----	4051642650241
	40351403	4051642005980	4051642150253	-----	4051642650258
	40351404	4051642005997	4051642150260	-----	4051642650265

The numerical data for the measurements, weights, thicknesses and tolerances are target values, and there may be some deviation.
The information represents the current state of our knowledge. It offers no guarantee of the product characteristics and does not form the basis of a contractual legal understanding. Please reconfirm any information that is older than 2 years.
Typographical and other errors excepted.
The product data sheets provided are the property of WIROS Wilfried Rosbach GmbH. WIROS Wilfried Rosbach GmbH is the owner of all copyrights and rights of use and exploitation for the product data sheets. It is not permitted to change the product data sheets or to use the pictograms. The product data sheets provided may only be used in conjunction with original WIROS products.

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